

## Protocol C-circle assay - AG Rippe

### 1. DNA concentration

- Measure DNA concentration with Qubit dsDNA HS Assay (Life Technologies)
- Dilute samples to a concentration of 10-30 ng/ $\mu$ l

### 2. Amplification

- Optional: Digest DNA with HinfI/RsaI and remove RNA by RNase digest
- Prepare 30 ng of each sample in 10  $\mu$ l H<sub>2</sub>O (every sample twice to include no Pol controls for sample)
- Prepare Master mix (without addition of polymerase first, see below)

	Stock	Final conc.	Vol per sample
BSA (NEB)	10 mg/ml	0.2 mg/ml	0.2 $\mu$ l
Tween 20	10%	0.1 %	0.1 $\mu$ l
dATP, dGTP, dTTP (Invitrogen)	100 mM	1 mM each	0.1 $\mu$ l each
$\Phi$ 29 buffer (NEB)	10x	2x	2 $\mu$ l
Nuclease-free water		Ad 10 $\mu$ l	6.65 $\mu$ l
$\Phi$ 29 polymerase (NEB)	10 U/ $\mu$ l	7.5 U	0.75 $\mu$ l

- **Before addition of polymerase: add 9.25  $\mu$ l Master mix to no Pol controls!**
- Add 10  $\mu$ l Master mix with polymerase to all other samples and mix
- In PCR thermocycler:
  - o 8 h @ 30°C
  - o 20 min @ 65°C
  - o 4°C hold

### 3. Dot Blot with Dot blotter (Bio Rad)

- Add 60  $\mu$ l 2X SSC to amplified samples in PCR tubes (total volume 80  $\mu$ l)
- Cut Nylon membrane (Roti-Nylon plus, 0.45  $\mu$ m) and soak 5 min in 2X SSC
- Place membrane on dot blot apparatus without trapping air bubbles!
- Load ~100  $\mu$ l 2X SSC in all wells and apply vacuum
- Load samples and fill empty wells with 2X SSC
- Apply vacuum
- Reload all wells with 100  $\mu$ l 2X SSC and apply vacuum again

- Let membrane dry for 10 min
- Bake membrane on whatman paper in glass dish for 20 min at 120 °C
- Wash membrane once in 2X SSC
- At this step the membrane can be dried and stored at 4 °C for several days
- Hybridization with DIG-telomere probe according to Telo TAGGG telomere length kit (Roche); in short (for each washing step prepare ~30 ml):
  - Prewarm DIG Easy Hyb Granules (bottle 7) to 42 °C
  - Prehybridization: Cover the membrane with prewarmed DIG Easy Hyb Granules in a tupper box and incubate for 30-60 min at 42 °C with gentle agitation
  - Prepare 5-10 ml hybridization solution depending on size of membrane (1 µl telomere probe per 5 ml prewarmed DIG Easy Hyb Granules)
  - Discard prehybridization solution and immediately add hybridization solution to the membrane
  - Incubate for 3 h (or O/N) at 42 °C with gentle agitation
  - Discard hybridization solution
  - Wash membrane twice for 5 min at RT with stringent wash buffer I (2X SSC, 0.1% SDS)
  - Wash membrane twice for 15-20 min at 50 °C with prewarmed stringent wash buffer II (0.2X SSC, 0.1% SDS) with gentle agitation
  - Rinse membrane in at least 30 ml 1X washing buffer for 5 min at RT with gentle agitation
  - Prepare ~60 ml 1X blocking solution: Dilute maleic acid buffer (bottle 12) 1:10 with H<sub>2</sub>O to get a 1X solution. Then, dilute 10X blocking buffer (bottle 11) 1:10 with 1X maleic acid buffer.
  - Incubate membrane in 30 ml freshly prepared 1X blocking solution for 30 min at RT with gentle agitation
  - Prepare Anti-DIG-AP working solution: Spin vial for 5 min at 13,000 rpm before use. Dilute Anti-DIG-AP (bottle 13) with 1X blocking solution to a final concentration of 75 mU/ml (1:10,000)
  - Incubate the membrane in 30 ml Anti-DIG-AP working solution for 30 min at RT with gentle agitation
  - Wash membrane 2x 15 min at RT with 1X washing buffer (bottle 10; 50 ml each time) with gentle agitation

- Prepare 1X detection buffer: Dilute 10X detection buffer (bottle 14) 1:10 with H<sub>2</sub>O
- Incubate membrane in 30 ml of 1X detection buffer for 2-5 min at RT with gentle agitation
- Discard detection buffer and remove excess liquid by placing the membrane on a whatman paper. Don't let the membrane dry!
- Place the wet membrane, DNA side facing up, on an open transparent sheet and very quickly apply about 2 ml substrate solution (bottle 15) to the membrane
- Immediately cover the membrane with the transparent sheet and spread the substrate solution without trapping air bubbles
- Incubate the membrane for 30 min to several hours at RT
- Squeeze out excess substrate solution
- Expose the membrane in a Bio-Rad Chemi-Doc MP System (Program "Chemi High Resolution")